



# Chromatography Vials

## Data sheet

### Crimp top vials and caps

- Large opening vials with 11mm crimp neck finish
- 1.8mL volume, standard size of 12x32mm diameter x height
- Excellent for long term storage
- Clear vials made from Type I Class A glass (33 expansion)
- Amber vials made from Type I Class B glass (for light sensitive product)
- Requires crimper (and decapper)-see accessories
- Sold in packs of 100 or 1000 as stated
- Manufactured in the USA and/or EC to precise standards
- 100% guaranteed



### Vials, 1.8mL, 12x32mm

Code	Description	List Price	Stock Measure
130297	Vial 1.8mL 12x32mm clear crimp	£76.96	1000pack
130298	Vial 1.8mL 12x32mm clear crimp w/spot	£62.07	1000pack
130299	Vial 1.8mL 12x32mm amber crimp	POA	1000pack
130300	Vial 1.8mL 12x32mm amber crimp w/spot	£100.00	1000pack
130297D	Vial 1.8mL 12x32mm clear crimp silanized	£21.60	100pack
130298D	Vial 1.8mL 12x32mm clr crimp w/spot sila	£22.93	100pack
130299D	Vial 1.8mL 12x32mm amber crimp silanized	£22.00	100pack
130300D	Vial 1.8mL 12x32mm amb crimp w/spot sila	£23.33	100pack

### Crimp Caps, 11mm, aluminium with septa

Code	Description	List Price	Stock Measure
130270	Cap crimp 11mm alum PTFE/red rub natural	£56.00	1000pack
130270B	Cap crimp 11mm alum PTFE/red rub blue	£82.66	1000pack
130270G	Cap crimp 11mm alum PTFE/red rub green	£82.66	1000pack
130270R	Cap crimp 11mm alum PTFE/red rub red	£82.66	1000pack
130275	Cap crimp 11mm alum PTFE/sil natural	£149.34	1000pack
130275B	Cap crimp 11mm alum PTFE/sil blue	£149.34	1000pack
130275G	Cap crimp 11mm alum PTFE/sil green	£149.34	1000pack
130275R	Cap crimp 11mm alum PTFE/sil red	£149.34	1000pack
130280	Cap crimp 11mm alum PTFE/sil/PTFE natur	£17.33	100pack
130280B	Cap crimp 11mm alum PTFE/sil/PTFE blue	£17.33	100pack
130280G	Cap crimp 11mm alum PTFE/sil/PTFE green	£17.33	100pack
130280R	Cap crimp 11mm alum PTFE/sil/PTFE red	£17.33	100pack



Distributed by:

Stratlab Ltd., Unit 1, Norwood Road, March, Cambridgeshire, PE15 8QD, United Kingdom

Tel.: +44 (0) 1354 657000 Fax.: +44 (0) 1354 656675

sales@stratlab.co.uk

www.stratlab.co.uk

For more information, visit: **www.dot-red.com**



## Chromatography Vials

### Data sheet

#### Snap/Crimp top vials

- Large opening vials with 11mm snap/crimp neck finish
- 1.8mL volume, standard size of 12x32mm diameter x height
- Clear vials made from Type I Class A glass (33 expansion)
- Amber vials made from Type I Class B glass (for light sensitive product)
- Accepts both crimp and snap caps (use crimper (see accessories) for crimp caps)
- Sold in packs of 100 or 1000 as stated
- See crimp top vials (11mm) for crimp caps
- Manufactured in the USA and/or EC to precise standards
- 100% guaranteed

#### Vials, 1.8mL, 12x32mm

Code	Description	List Price	Stock Measure
130291	Vial 1.8mL 12x32mm clear snap	£80.87	1000pack
130292	Vial 1.8mL 12x32mm clear snap w/spot	£89.13	1000pack
130293	Vial 1.8mL 12x32mm amber snap	POA	1000pack
130294	Vial 1.8mL 12x32mm amber snap w/spot	£93.92	1000pack
130291D	Vial 1.8mL 12x32mm clear snap silanized	£21.60	100pack
130292D	Vial 1.8mL 12x32mm clr snap w/spot silan	£22.93	100pack
130293D	Vial 1.8mL 12x32mm amber snap silanized	£22.00	100pack
130294D	Vial 1.8mL 12x32mm amb snap w/spot silan	£23.33	100pack



Distributed by:

Stratlab Ltd., Unit 1, Norwood Road, March, Cambridgeshire, PE15 8QD, United Kingdom

Tel.: +44 (0) 1354 657000 Fax.: +44 (0) 1354 656675


sales@stratlab.co.uk

www.stratlab.co.uk

For more information, visit: **www.dot-red.com**

### Snap Caps, 11mm, polyethylene, with septa

Code	Description	List Price	Stock Measure
130323	Cap snap 11mm PE PTFE/red rubber natural	£122.66	1000pack
130323B	Cap snap 11mm PE PTFE/red rubber blue	£120.00	1000pack
130323BK	Cap snap 11mm PE PTFE/red rubber black	£122.66	1000pack
130323G	Cap snap 11mm PE PTFE/red rubber green	£122.66	1000pack
130323R	Cap snap 11mm PE PTFE/red rubber red	£122.66	1000pack
130323Y	Cap snap 11mm PE PTFE/red rubber yellow	£122.66	1000pack
130328	Cap snap 11mm PE PTFE sil natural	£161.34	1000pack
130328B	Cap snap 11mm PE PTFE sil blue	£161.34	1000pack
130328BK	Cap snap 11mm PE PTFE sil black	£161.34	1000pack
130328G	Cap snap 11mm PE PTFE sil green	£161.34	1000pack
130328R	Cap snap 11mm PE PTFE sil red	£161.34	1000pack
130328Y	Cap snap 11mm PE PTFE sil yellow	£161.34	1000pack
130333	Cap snap 11mm PE PTFE/sil/PTFE natural	£20.40	100pack
130333B	Cap snap 11mm PE PTFE/sil/PTFE blue	£20.40	100pack
130333BK	Cap snap 11mm PE PTFE/sil/PTFE black	£20.40	100pack
130333G	Cap snap 11mm PE PTFE/sil/PTFE green	£20.40	100pack
130333R	Cap snap 11mm PE PTFE/sil/PTFE red	£20.40	100pack
130333Y	Cap snap 11mm PE PTFE/sil/PTFE yellow	£20.40	100pack
130334	Cap snap 11mm PE PTFE/sil w/slit natural	£23.33	100pack
130334B	Cap snap 11mm PE PTFE/sil w/slit blue	£23.33	100pack
130334BK	Cap snap 11mm PE PTFE/sil w/slit blk	£23.33	100pack
130334G	Cap snap 11mm PE PTFE/sil w/slit green	£23.33	100pack
130334R	Cap snap 11mm PE PTFE/sil w/slit red	£23.33	100pack
130334Y	Cap snap 11mm PE PTFE/sil w/slit yellow	£23.33	100pack



Distributed by:

Stratlab Ltd., Unit 1, Norwood Road, March, Cambridgeshire, PE15 8QD, United Kingdom

Tel.: +44 (0) 1354 657000 Fax.: +44 (0) 1354 656675

sales@stratlab.co.uk

www.stratlab.co.uk

For more information, visit: **www.dot-red.com**



## Chromatography Vials

### Data sheet

#### Screw thread vials

- Large neck opening vials with 9mm neck finish
- 1.8mL volume, standard size of 12x32mm diameter x height
- Excellent for long term sample storage
- Clear vials made from Type I Class A glass (33 expansion)
- Amber vials made from Type I Class B glass (for light sensitive product)
- Sold in packs of 100 or 1000 as stated
- Manufactured to precise standards in the USA and/or EC
- 100% guaranteed

#### Vials, 1.8mL, 12 x 32mm

Code	Description	List Price	Stock Measure
130638	Vial 1.8mL 12x32mm clear screw 9mm	£101.34	1000pack
130639	Vial 1.8mL 12x32mm clr screw w/spot 9mm	£111.10	1000pack
130640	Vial 1.8mL 12x32mm amber screw 9mm	POA	1000pack
130641	Vial 1.8mL 12x32mm amber scrw w/spot 9mm	£108.00	1000pack
130638D	Vial 1.8mL 12x32mm clear screw silan 9mm	£22.27	100pack
130639D	Vial 1.8mL 12x32mm clr scrw+spot sil 9mm	£23.60	100pack
130640D	Vial 1.8mL 12x32mm amber scrw silan 9mm	£22.80	100pack
130641D	Vial 1.8mL 12x32mm amb scrw+spot sil 9mm	£24.13	100pack



Distributed by:

Stratlab Ltd., Unit 1, Norwood Road, March, Cambridgeshire, PE15 8QD, United Kingdom

Tel.: +44 (0) 1354 657000 Fax.: +44 (0) 1354 656675


sales@stratlab.co.uk

www.stratlab.co.uk

For more information, visit: **www.dot-red.com**

### Screw thread caps, 9mm, polypropylene, with septa

Code	Description	List Price	Stock Measure
130241B	Cap screw 9mm PE PTFE/red rubber blue	£110.87	1000pack
130241BK	Cap screw 9mm PE PTFE/red rubber black	£132.00	1000pack
130241G	Cap screw 9mm PE PTFE/red rubber green	£132.00	1000pack
130241R	Cap screw 9mm PE PTFE/red rubber red	£132.00	1000pack
130241Y	Cap screw 9mm PE PTFE/red rubber yellow	£132.00	1000pack
130242B	Cap screw 9mm PE PTFE/silicone blue	£146.66	1000pack
130242BK	Cap screw 9mm PE PTFE/silicone black	£129.13	1000pack
130242G	Cap screw 9mm PE PTFE/silicone green	£146.66	1000pack
130242R	Cap screw 9mm PE PTFE/silicone red	£146.66	1000pack
130242Y	Cap screw 9mm PE PTFE/silicone yellow	£146.66	1000pack
130243B	Cap screw 9mm PE PTFE/silicon/PTFE blue	£18.00	100pack
130243BK	Cap screw 9mm PE PTFE/silicon/PTFE black	£18.00	100pack
130243G	Cap screw 9mm PE PTFE/silicon/PTFE green	£18.00	100pack
130243R	Cap screw 9mm PE PTFE/silicon/PTFE red	£18.00	100pack
130243Y	Cap screw 9mm PE PTFE/silic/PTFE yellow	£18.00	100pack
130244B	Cap screw 9mm PE PTFE/sili w/slit blue	£22.53	100pack
130244BK	Cap screw 9mm PE PTFE/sili w/slit black	£22.53	100pack
130244G	Cap screw 9mm PE PTFE/sil w/slit green	£22.53	100pack
130244R	Cap screw 9mm PE PTFE/silicon w/slit red	£22.53	100pack
130244Y	Cap screw 9mm PE PTFE/sil w/slit yellow	£22.53	100pack



Distributed by:

Stratlab Ltd., Unit 1, Norwood Road, March, Cambridgeshire, PE15 8QD, United Kingdom

Tel.: +44 (0) 1354 657000 Fax.: +44 (0) 1354 656675

sales@stratlab.co.uk

www.stratlab.co.uk

For more information, visit: **www.dot-red.com**

### Head space vials and caps

- 20mm crimp cap neck finish
- Clear vials manufactured in Type I Class A glass (33 Expansion)
- Pressure release seals offer a safety feature where the scored aluminium seal will open if the internal pressure reaches  $3.0 \pm 0.5$  bar (1 bar = ~ 14.5 psi)
- Sold in packs of 100
- Manufactured in the USA and/or EC to precise standards
- 100% guaranteed



### Vials, 6, 10 and 20mL

Code	Description	List Price	Stock Measure
130301	Vial 6mL 22x38mm clear flat bot std neck	£16.13	100pack
130302	Vial 10mL 23x46mm clear flat bot std nec	£18.43	100pack
130303	Vial 20mL 23x75mm clear flat bot std nec	£19.33	100pack
130304	Vial 6mL 22x38mm clr flat bot taper neck	£15.78	100pack
130305	Vial 10mL 23x46mm clr flat bot taper nec	£18.27	100pack
130306	Vial 20mL 23x75mm clr flat bot taper nec	£19.33	100pack
130307	Vial 6mL 23x38mm clr rd bot taper neck	£16.13	100pack
130308	Vial 10mL 23x46mm clr rd bot taper neck	£18.27	100pack
130309	Vial 20mL 23x75mm clr rd bot taper neck	£19.33	100pack

### Crimp caps, 20mm aluminium (unless stated) with septa

Code	Description	List Price	Stock Measure
130315	Cap crimp 20mm alum PTFE/silicone	£27.22	100pack
130316	Cap crimp 20mm alum PTFE/butyl rubber	£37.20	100pack
130317	Cap crimp 20mm alum PTFE/moulded butyl	£4.24	100pack
130318	Cap crimp 20mm alum PTFE/sil pres rel	£42.00	100pack
130319	Cap crimp 20mm alum PTFE/butyl pres rel	£50.93	100pack
130320	Cap crimp 20mm al PTFE/mouldbutyl presre	£55.87	100pack
130321	Cap crimp 20mm PTFE/silicone magnetic	£44.93	100pack
130322	Cap crimp 20mm PTFE/butyl magnetic	£55.60	100pack
130326	Cap crimp 20mm PTFE/mouldbutyl magnetic	£57.73	100pack

Distributed by:

Stratlab Ltd., Unit 1, Norwood Road, March, Cambridgeshire, PE15 8QD, United Kingdom

Tel.: +44 (0) 1354 657000 Fax.: +44 (0) 1354 656675

sales@stratlab.co.uk

www.stratlab.co.uk

For more information, visit: **www.dot-red.com**

### Accessories

#### Crimpers and Decappers

Code	Description	List Price	Stock Measure
130475	Crimper hand 11mm for 12mm cap vial	£117.99	each
130476	Decapper hand op crimper type 11mm cap	£134.34	each
130528	Decapper hand op plier type 11mm cap	£34.12	each
130587	Crimper hand 20mm for 20mm cap vial	£117.99	each
130587F	Crimper hand 20mm Ø 23mm cap vial flip	£118.30	each
130589	Decapper hand op crimper type 20mm cap	£132.34	each
130592	Decapper hand op plier type 20mm cap	£33.52	each



#### Inserts, 0.25mL

Code	Description	List Price	Stock Measure
130324	Insert 0.25mL glass w/bottom spr 6x29mm	£386.66	100pack
130325	Insert 0.25mL flat bottom 6x31mm	£50.00	1000pack

Distributed by:

Stratlab Ltd., Unit 1, Norwood Road, March, Cambridgeshire, PE15 8QD, United Kingdom

Tel.: +44 (0) 1354 657000 Fax.: +44 (0) 1354 656675

sales@stratlab.co.uk

www.stratlab.co.uk

For more information, visit: **www.dot-red.com**



## Chromatography Vials

### Data sheet

#### Other vials

#### Combinatorial chemistry vials

- Available in a variety of sizes and styles
- Custom vial designs available to fit your automation requirements
- Custom bar coding services-ceramic and paper
- Closures available for all sizes
- Custom engineered rack designs available
- Manufactured in the USA and/or EC to precise standards
- 100% guaranteed

Code	Description	List Price	Stock Measure
130312	Vial 0.3mL 12x32mm clear crimp comb 11mm	POA	100pack
130313	Vial 0.3mL 12x32mm clr snap comb 11mm	POA	100pack
130314	Vial 0.3mL 12x32mm clear screw comb 9mm	POA	100pack
130310	Vial 3.1mL 15x45mm clear screw comb	£138.26	100pack
130311	Vial 3.1mL 15x51mm clear screw comb	POA	100pack
130329	Vial 5mL 17x60mm clear screw comb	POA	100pack
130330	Vial 10mL 21x70mm clear screw comb	POA	100pack
130331	Vial 20mL 28x67mm clear screw comb	POA	100pack
130332	Vial 33mL 28x88mm clear screw comb	POA	100pack

Distributed by:

Stratlab Ltd., Unit 1, Norwood Road, March, Cambridgeshire, PE15 8QD, United Kingdom

Tel.: +44 (0) 1354 657000 Fax.: +44 (0) 1354 656675

sales@stratlab.co.uk

www.stratlab.co.uk

For more information, visit: **www.dot-red.com**